



AA200892

Long-life recirculation kit, initial equipment

AA200892

- For air recirculation.
- Easy removal of the air recirculation module from the top of the table ventilation for regenerating.
- Regenerable in the oven at 200 °C.
- Extended change interval of the activated charcoal filter: regenerable up to 30 times.
- For air recirculation.

Easy removal of the air recirculation module from the top of the table ventilation for regenerating.

Regenerable in the oven at 200 °C.

Extended change interval of the activated charcoal filter: regenerable up to 30 times.

